

Product Datasheet

Anti-Atrogin-1 Antibody (orb1822312)

Catalog Number	orb1822312
Category	Antibodies
Description	Our Atrogin-1 rabbit polyclonal primary antibody from Biorbyt is produced in-house. It detects human, mouse, and rat Atrogin-1 and is antigen affinity purified.
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Form/Appearance	Liquid
Buffer/Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Immunogen	A synthetic peptide (coupled to carrier protein) corresponding to amino acids from the N-terminal region of mouse Atrogin. This sequence is highly homologous to human and rat Atrogin (F-box only protein 32).
UniProt ID	Q9CPU7
MW	41
Tested applications	IHC, WB
Dilution range	WB: 1:1000, IHC: 1:100
Antibody Type	Primary Antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)