
Product Datasheet

Anti-Arpc1b/p41-Arc (Thr-21), Phosphospecific Antibody (orb1822318)

Description

Anti-Arpc1b/p41-Arc (Thr-21),

Species/Host	Rabbit
Reactivity	Gallus, Human, Mouse, Rat, Xenopus
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Phospho-Arpc1b/p41-Arc (Thr-21) synthetic peptide (coupled to KLH) corresponding to amino acid residues surrounding threonine 21 in human Arpc1b. This peptide sequence is highly conserved in rat, mouse, chicken, and fish Arpc1b, and has some homology to the conserved site in isoform Arpc1a.
Preservatives	PBS + 1 mg/ml rAlbumin, 0.05% NaN ₃ and 50% glycerol
Form/Appearance	Antigen Affinity Purified
Storage	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41
Uniprot ID	015143
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.