

---

## Product Datasheet

### Anti-Arpc1b/p41-Arc (C-terminal region) Antibody (orb1822319)

**Description**

Anti-Arpc1b/p41-Arc (C-terminal)

**Species/Host**

Rabbit

**Reactivity**

Gallus, Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

ICC, IP, WB

**Immunogen**

Arpc1b/p41-Arc synthetic peptide (coupled to KLH) corresponding to amino acid residues in the C-terminal region of human Arpc1b. This peptide sequence is highly conserved in rat, mouse, and chicken Arpc1b, and has some homology to the conserved region in isoform Arpc1a.

**Preservatives**PBS + 1 mg/ml rAlbumin, 0.05% NaN<sub>3</sub> and 50% glycerol**Form/Appearance**

Antigen Affinity Purified

**Storage**

Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**MW**

41

**Uniprot ID**[O15143](#)**Dilution Range**

WB: 1:1000, ICC: 1:200, IP: 1:100

**Expiration Date**

12 months from date of receipt.