

---

## Product Datasheet

### HSP90 beta Protein (orb1822409)

**Description**

Active Human Recombinant HSP90 beta Full Length Protein

**Reactivity**

Human

**Tested Applications**

CoIP, ELISA, FA, SBR, SDS-PAGE, WB

**Target**

HSP90 beta

**Preservatives**

20mM Tris, pH 7.5, 175 mM NaCl, 0.1 mM EDTA, 10% glycerol, 1 mM DTT

**Concentration**

Lot/batch specific. See included datasheet.

**Storage**

-20°C

**Tag**

No tag

**Note**

For research use only

**Application notes**

This product has been certified &gt; 90% pure using SDS-PAGE analysis.

**Protein Sequence**

MPEEVHHGEEEVETFAFQAEIAQLMSLIINTFYNSKEIFLRELISNASDALDKIRYES  
 LTDPSKLDGKELKIDIIPNPQERTLTLVDTGIGMTKADLNNLGTIAKSGTKAFME  
 ALQAGADISMIGQFGVGFYSAYLVAEKVVITKHNDDEQYAWESSAGGSFTVRA  
 DHGEPGRGTVILHLKEDQTEYLEERRVKEVVKHSQFIGYPITLYLEKEREKEIS  
 DDEAEKEEKEEEDKDEEKPKIEDVGSDEEDDSGKDKKKKTKKIKEKYIDQEE  
 LNKTPIWTRNPDDITQEEYGEFYKSLTNDWEDHLAVKHFSVEGQLEFRALLFIP  
 RRAPFDLFENKKKNNIKLYVRRVFMDSCELIPEYLNFIKRVVDSDELPLNISRE  
 MLQQSKILKVIKKNIVKKCLELFSCLAEDKENYKFFYAFSKNLKLGIEDSTNRRR  
 LSELLRYHTSQSGDEMTSLSEYVSRMKETQKSIYYITGESKEQVANSFVERVRK  
 RGFVYVMTPEIDEYCVQLKEFDGKSLVSVTKEGLELPEDEEKKKMEESKAKF  
 ENLCKLMKEILDKKVEKVTISNRLVSSPCCIVTSTYGWTANMERIMKAQALRDN  
 TMGYMMAKKHLEINPDHPVETLRQKAEADKNDKAVKDLVLLFETALLSSGFSL  
 EDPQTHSNRIYRMIKGLGLIDEDEVAEEPNAAVPDEIPPLEGDEDEDASRMEEVD

**Purity**

&gt; 85%

**Source**

Recombinant

**MW**

~90 kDa

**Uniprot ID**
**P08238**
**Expiration Date**

6 months from date of receipt.