



Product Datasheet

ATRX Antibody / RAD54 (orb1822538)

Description	ATRX is a member of the Snf2 family of helicase/ATPases, which contribute to the remodeling of the n...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human ATRX protein was used as the immunogen for the ATRX antibody.
Storage	Aliquot the ATRX antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	P46100
Hazard Information	This ATRX antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon tissue ...



SDS-PAGE analysis of purified, BSA-free ...