



Product Datasheet

Androgen Receptor Antibody (orb1822542)

Description	Androgen Receptor is a member of the superfamily of ligand responsive transcription regulators. The ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A synthetic peptide corresponding to a portion of human Androgen Receptor protein (within amino acids 1-100) was used as the immunogen for the Androgen Receptor antibody.
Storage	Aliquot the Androgen Receptor antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P10275
Hazard Information	This Androgen Receptor antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate carc...



IHC staining of FFPE human breast carino...