



Product Datasheet

ACE2 / Angiotensin I Converting Enzyme 2 Antibody (orb1822595)

Description	Carboxypeptidase which converts angiotensin I to angiotensin 1-9, a peptide of unknown function, and...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-300) from the human protein was used as the immunogen for the ACE2 antibody.
Storage	Aliquot the ACE2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q9BYF1
Hazard Information	This ACE2 antibody is available for research use only.
Dilution Range	ELISA (Order BSA-free format for coating), Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human small intesti...



IHC staining of FFPE human kidney tissue...