

Product Datasheet

Anti-HSP60 Hspd1 Antibody (Monoclonal, LK1) (orb18226)

Description	Anti-HSP60 Hspd1 Antibody (Monoclonal, LK1)
Species/Host	Mouse
Reactivity	Gallus, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant human heat shock protein 60 (HSP60)
Form/Appearance	Lyophilized
Concentration	Adding 1 ml of PBS buffer will yield a concentration of 100 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Western blot, 2-4µg/ml, Human, chicken, rat Immunohistochemistry (Paraffin-embedded Section), 4-8µg/ml, Human, chicken, mouse, rat. Add 1ml of PBS buffer will yield a concentration of 100ug/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Clone Number	LK1
MW	65 kDa
Uniprot ID	P63039


Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Expiration Date 12 months from date of receipt.



WB analysis of HSP60 using anti-HSP60 antibody. Lane 1: human HeLa cell; 2: human U87 cell; 3: rat brain tissue; 4: mouse Neuro-2a cell.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)