



Product Datasheet

TFE3 Antibody / Transcription factor E3 (orb1822637)

Description	Transcription factor E3 (TFE3) belongs to the basic helix-loop-helix zipper transcription factor fam...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human TFE3 protein was used as the immunogen for the TFE3 antibody.
Storage	Aliquot the TFE3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P19532
Hazard Information	This TFE3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human MIT transloca...



Negative control: IHC testing of FFPE hu...