


Product Datasheet

TFE3 Antibody / Transcription factor E3 (orb1822640)

Description	Transcription factor E3 (TFE3) belongs to the basic helix-loop-helix zipper transcription factor fam...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human protein was used as the immunogen for the Transcription factor E3 antibody.
Storage	Aliquot the Transcription factor E3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P19532
Hazard Information	This Transcription factor E3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thyroid tissue...