


Product Datasheet

TFE3 Antibody / Transcription factor E3 (orb1822641)

Description	Transcription factor E3 (TFE3) belongs to the basic helix-loop-helix zipper transcription factor fam...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human protein was used as the immunogen for the Transcription factor E3 antibody.
Storage	Aliquot the Transcription factor E3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P19532
Hazard Information	This Transcription factor E3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human thyroid tissu...