

Product Datasheet

KDM5C Antibody / JARID1C / SMCX (orb1822644)

Description SmcX, also known as JARID1C (jumonji, AT rich interactive domain 1C), MRXJ, KDM5C or XE169, is a nuc...

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FACS

Immunogen Recombinant full-length human protein was used as the immunogen for the JARID1C antibody.

Storage Aliquot the JARID1C antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note For research use only

Formula 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype Mouse IgG2c


Clonality Monoclonal

Uniprot ID **P41229**


Hazard Information This JARID1C antibody is available for research use only.

Dilution Range Flow cytometry: 1-2ug/million cells


Expiration Date 12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



Flow cytometry testing of PFA-fixed huma...



SDS-PAGE analysis of purified, BSA-free ...