


## Product Datasheet

# HSV1 Antibody / Herpes Simplex Virus Type I (orb1822647)

<b>Description</b>	The antibody reacts with HSV type 1 specific antigen. It is suitable for detection of HSV in human C...
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	Detergent-solubilized herpes simplex virus (HSV) type 1 infected cells were used as the immunogen for the Herpes Simplex Virus Type I antibody.
<b>Storage</b>	Aliquot the Herpes Simplex Virus Type I antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
<b>Isotype</b>	Rabbit IgG, kappa
<b>Clonality</b>	Recombinant
<b>Hazard Information</b>	This Herpes Simplex Virus Type I antibody is available for research use only.
<b>Dilution Range</b>	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
<b>Expiration Date</b>	12 months from date of receipt.



IHC staining of FFPE human HSV-infected ...