



Product Datasheet

Helicobacter pylori Catalase Antibody (orb1822654)

Description	The spiral shaped bacterium Helicobacter pylori is strongly associated with inflammation of the stom...
Species/Host	Rabbit
Reactivity	H. pylori
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant Helicobacter pylori Catalase protein fragment (within amino acids 323-445) was used as the immunogen for the Helicobacter pylori antibody.
Storage	Aliquot the Helicobacter pylori antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Hazard Information	This Helicobacter pylori antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human H. pylori-inf...



SDS-PAGE analysis of purified, BSA-free ...