



Product Datasheet

HBsAg Antibody / Hepatitis B surface antigen (orb1822669)

Description	Hepatitis B surface antigen (HBsAg) is a glycoprotein on the surface of hepatitis B virus. After hep...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human HBsAg protein was used as the immunogen for the HBsAg antibody.
Storage	Aliquot the HBsAg antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P03138
Hazard Information	This HBsAg antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human HBV-infected ...



SDS-PAGE analysis of purified, BSA-free ...