

Product Datasheet

HBsAg Antibody / Hepatitis B surface antigen (orb1822670)

Description	Hepatitis B surface antigen (HBsAg) is a glycoprotein on the surface of hepatitis B virus. After hepa...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human HBsAg protein was used as the immunogen for the Hepatitis B surface antigen antibody.
Storage	Aliquot the Hepatitis B surface antigen antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P03138
Hazard Information	This Hepatitis B surface antigen antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2 µg/ml for 30 min at RT
Expiration Date	12 months from date of receipt.