



Product Datasheet

CD105 Antibody / Endoglin (orb1822694)

Description	CD105/Endoglin is a Type I membrane glycoprotein located on cell surfaces and is part of the TGF-beta...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant extracellular fragment of human Endoglin protein was used as the immunogen for the Endoglin antibody.
Storage	Aliquot the Endoglin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P17813
Hazard Information	This Endoglin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human placental tis...



SDS-PAGE analysis of purified, BSA-free ...