
Product Datasheet

**Steroidogenic Factor 1 Antibody / SF-1 / NR5A1
(orb1822697)**

Description

Steroidogenic Factor 1 (SF-1) is considered an

5. ...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

A recombinant partial protein sequence (within amino acids 200-400) from the human protein was used as the immunogen for the SF-1 antibody.

Storage

Aliquot the SF-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Uniprot ID

[Q13285](#)

Hazard Information


This SF-1 antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



IHC staining of FFPE human testis tissue...



IHC staining of FFPE human testis tissue...



SDS-PAGE analysis of purified, BSA-free ...