



Product Datasheet

Steroidogenic Factor 1 Antibody / SF-1 / NR5A1 (orb1822702)

Description	Steroidogenic Factor 1 (SF-1) is considered an orphan nuclear receptor that belongs to subfamily 5. ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 200-400) from the human protein was used as the immunogen for the Steroidogenic Factor 1 antibody.
Storage	Aliquot the Steroidogenic Factor 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q13285
Hazard Information	This Steroidogenic Factor 1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human adrenal corte...



IHC staining of FFPE human testis tissue...