



Product Datasheet

TNFSF15 Antibody / VEGI / TL1A (orb1822712)

Description	VEGI is an anti-angiogenic cytokine that belongs to tumor necrosis factor superfamily, member 15 (TN...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human protein was used as the immunogen for the TNFSF15 antibody.
Storage	Aliquot the TNFSF15 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	O95150
Hazard Information	This TNFSF15 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human colon tissue
...



SDS-PAGE analysis of purified, BSA-free ...