



Product Datasheet

Alpha 1 microglobulin Antibody / AMBP / Bikunin (orb1822715)

Description	The AMBP (Alpha-1-Microglobulin/Bikunin precursor) gene encodes a protein precursor, known as AMBP, ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 200-353) from the human protein was used as the immunogen for the Alpha 1 Microglobulin antibody.
Storage	Aliquot the Alpha 1 Microglobulin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	P02760
Hazard Information	This Alpha 1 Microglobulin antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human placental tis...



SDS-PAGE analysis of purified, BSA-free ...