
Product Datasheet

**Alpha 1 microglobulin Antibody / AMBP / Bikunin
(orb1822718)**

Description

The AMBP (Alpha-1-Microglobulin/Bikunin

known as AMBP, ...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

A recombinant partial protein sequence (within amino acids 200-353) from the human protein was used as the immunogen for the AMBP antibody.

Storage

Aliquot the AMBP antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Mouse IgG2, kappa

Clonality

Monoclonal

Uniprot ID
P02760
Hazard Information


This AMBP antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



SDS-PAGE analysis of purified, BSA-free ...



IHC staining of FFPE human ovarian cancer...