



Product Datasheet

TLE1 Antibody (orb1822738)

Description	Key players in the Notch pathway are the TLE genes, which are human homologs of the Drosophila grouc...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant human protein fragment (aa 175-338) was used as the immunogen for the Transducin-like enhancer protein 1 antibody.
Storage	Aliquot the Transducin-like enhancer protein 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q04724
Hazard Information	This TLE1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human prostate carc...



IHC staining of FFPE human bladder tissu...