



Product Datasheet

TRP1 / Tyrosinase-Related Protein-1 Antibody (orb1822772)

Description	It reacts with a 75kDa melanocyte-specific gene product, identified as Tyrosinase-related protein-1 ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 200-400) from the human protein was used as the immunogen for the Tyrosinase-Related Protein-1 antibody.
Storage	Aliquot the Tyrosinase-Related Protein-1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P17643
Hazard Information	This Tyrosinase-Related Protein-1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human melanoma tiss...



SDS-PAGE analysis of purified, BSA-free ...