
Product Datasheet

**CA9 Antibody / Carbonic Anhydrase IX
(orb1822796)**

Description

Carbonic Anhydrase IX catalyzes the interconversion

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

A recombinant human protein fragment (within amino acids 314-410) was used as the immunogen for the Carbonic Anhydrase IX antibody.

Storage

Aliquot the Carbonic Anhydrase IX antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Uniprot ID

[Q16790](#)

Hazard Information


This Carbonic Anhydrase IX antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



IHC staining of FFPE human kidney tissue...



IHC staining of FFPE human kidney tissue...



IHC staining of FFPE human parathyroid t...