
Product Datasheet

**CA9 Antibody / Carbonic Anhydrase IX
(orb1822798)**

Description

Carbonic Anhydrase IX catalyzes the interconversion


Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

IHC-P




IHC staining of FFPE human kidney tissue...

Immunogen

A recombinant human protein fragment (within amino acids 314-410) was used as the immunogen for the Carbonic Anhydrase IX antibody.

Storage

Aliquot the Carbonic Anhydrase IX antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.



IHC staining of FFPE human kidney tissue...

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Uniprot ID


Q16790

Hazard Information

This Carbonic Anhydrase IX antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT



IHC staining of FFPE human parathyroid t...

Expiration Date

12 months from date of receipt.