

---

## Product Datasheet

**TRPS1 Antibody / Tricho-rhino-phalangeal  
syndrome type I (orb1822823)**

**Description**

Transcriptional repressor. Binds specifically to


**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P


IHC staining of FFPE human breast tissue...

**Immunogen**

A recombinant partial protein sequence (within amino acids 900-1100) from the human protein was used as the immunogen for the TRPS1 antibody.

**Storage**

Aliquot the TRPS1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.


**Note**

For research use only

**Formula**

1 mg/ml in 1X PBS; BSA free, sodium azide free

**Isotype**

Rabbit IgG, kappa

IHC staining of FFPE human mammary gland...

**Clonality**

Recombinant

**Uniprot ID**


**Q9UHF7**

**Hazard Information**

This TRPS1 antibody is available for research use only.

**Dilution Range**

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT


**Expiration Date**

12 months from date of receipt.

IHC staining of FFPE human uterus. Stron...