

---

## Product Datasheet

**TRPS1 Antibody / Tricho-rhino-phalangeal  
syndrome type I (orb1822824)**

**Description**

Transcriptional repressor. Binds specifically to


**Species/Host** Rabbit

**Reactivity** Human

**Conjugation** Unconjugated

**Tested Applications** IHC-P

**Immunogen** A recombinant partial protein sequence (within amino acids 900-1100) from the human protein was used as the immunogen for the TRPS1 antibody.

**Storage** Aliquot the TRPS1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

**Note** For research use only


**Formula** 0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

**Isotype** Rabbit IgG, kappa


**Clonality** Recombinant

**Uniprot ID** [Q9UHF7](#)
**Hazard Information** This TRPS1 antibody is available for research use only.


**Dilution Range** Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date** 12 months from date of receipt.


IHC staining of FFPE human breast tissue...



IHC staining of FFPE human mammary gland...



IHC staining of FFPE human uterus. Stron...