



Product Datasheet

TRPS1 Antibody / Tricho-rhino-phalangeal syndrome type I (orb1822831)

Description	Transcriptional repressor. Binds specifically to GATA sequences and represses expression of GATA-reg...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 900-1100) from the human protein was used as the immunogen for the Tricho-rhino-phalangeal syndrome type I antibody.
Storage	Aliquot the Tricho-rhino-phalangeal syndrome type I antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q9UHF7
Hazard Information	This Tricho-rhino-phalangeal syndrome type I antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human ovary tissue ...



IHC staining of FFPE human breast carcin...