



Product Datasheet

CA8 Antibody / Carbonic Anhydrase VIII (orb1822886)

Description	The protein encoded by this gene was initially named CA-related protein because of sequence similari...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human protein was used as the immunogen for the Carbonic Anhydrase VIII antibody.
Storage	Aliquot the Carbonic Anhydrase VIII antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	P35219
Hazard Information	This Carbonic Anhydrase VIII antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



IHC staining of FFPE human brain tissue

...