



Product Datasheet

GATA4 Antibody (orb1822889)


Description	Members of the GATA family share a conserved zinc finger DNA-binding domain and are capable of bindi...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FACS
Immunogen	A recombinant partial protein sequence (within amino acids 262-321) from the human protein was used as the immunogen for the GATA4 antibody.
Storage	Aliquot the GATA4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2b
Clonality	Monoclonal
Uniprot ID	P43694
Hazard Information	This GATA4 antibody is available for research use only.
Dilution Range	Flow cytometry: 1-2ug/million cells
Expiration Date	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



Flow cytometry testing of PFA-fixed huma...



SDS-PAGE analysis of purified, BSA-free ...