
Product Datasheet

**DBC1 Antibody / Deleted in breast cancer 1 /
CCAR2 (orb1822896)**

Description

DBC-1 (deleted in breast cancer gene 1 protein),

ap...

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FACS, IHC-P

Immunogen

Recombinant full-length human CCAR2 / KIAA1967 protein was used as the immunogen for the DBC1 antibody.

Storage

Aliquot the DBC1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a

Clonality

Monoclonal

Uniprot ID

Q8N163

Hazard Information


This DBC1 antibody is available for research use only.

Dilution Range


Flow cytometry: 1-2ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



Flow cytometry testing of PFA-fixed huma...



IHC staining of FFPE human colon tumor t...



Analysis of a HuProt (TM) microarray con...