



Product Datasheet

EZH2 Antibody (orb1822913)


Description	This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric pr...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A recombinant partial protein sequence (within amino acids 300-550) from the human protein was used as the immunogen for the EZH2 antibody.
Storage	Aliquot the EZH2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q15910
Hazard Information	This EZH2 antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human renal cell ca...



Western blot testing of human HDLM-2 cel...



SDS-PAGE analysis of purified, BSA-free ...