



Product Datasheet

EZH2 Antibody (orb1822918)


Description	This gene encodes a member of the Polycomb-group (PcG) family. PcG family members form multimeric pr...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A recombinant partial protein sequence (within amino acids 300-550) from the human protein was used as the immunogen for the Enhancer of zeste homolog 2 antibody.
Storage	Aliquot the EZH2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG, kappa
Clonality	Monoclonal
Uniprot ID	Q15910
Hazard Information	This Enhancer of zeste homolog 2 antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human bladder carcinoma tissue.



IHC staining of FFPE human Lynch syndrome tissue.



IHC staining of FFPE human bladder carcinoma tissue.