



Product Datasheet

EPHB4 Antibody / EPH Receptor B4 (orb1822979)


Description	The Eph subfamily represents the largest group of receptor protein tyrosine kinases identified to da...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the Ephrin type-B receptor 4 antibody.
Storage	Aliquot the Ephrin type-B receptor 4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	P54760
Hazard Information	This Ephrin type-B receptor 4 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human liver tissue ...



IHC staining of FFPE human salivary gland...



Analysis of a HuProt (TM) microarray con...