





Product Datasheet

EPHB4 Antibody / EPH Receptor B4 (orb1822986)

Description	The Eph subfamily represents the largest group of receptor protein tyrosine kinases identified to da...	
Species/Host	Mouse	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P	
Immunogen	A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the EPH Receptor B4 antibody.	
Storage	Aliquot the EPH Receptor B4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.	
Note	For research use only	
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	
Isotype	Mouse IgG2, kappa	
Clonality	Monoclonal	
Uniprot ID	P54760	
Hazard Information	This EPH Receptor B4 antibody is available for research use only.	
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT	
Expiration Date	12 months from date of receipt.	