



Product Datasheet

PON1 Antibody (orb1822988)


Description	Paroxon is an organophosphorus anticholinesterase compound, used topically in the treatment of glauc...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the Paraoxonase 1 antibody.
Storage	Aliquot the Paraoxonase 1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	P27169
Hazard Information	This Paraoxonase 1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human liver tissue ...



Analysis of a HuProt (TM) microarray con...



SDS-PAGE analysis of purified, BSA-free ...