


Product Datasheet

HSP27 Antibody / HSPB1 (orb1822997)

Description	HSP27, also referred to as the Estrogen-Regulated 24K protein and HSP28, is one of several small hea...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant human full-length protein was used as the immunogen for the HSPB1 antibody.
Storage	Aliquot the HSPB1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	1 mg/ml in 1X PBS; BSA free, sodium azide free
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P04792
Hazard Information	This HSPB1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human renal cell ca...