



Product Datasheet

IL22RA2 Antibody (orb1823061)

Description	IL-22RA α 2 (Interleukin-22 receptor Alpha2 chain), also known as IL-22 binding protein (IL-22BP), ...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-200) from the human protein was used as the immunogen for the IL22RA2 antibody.
Storage	Aliquot the IL22RA2 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2, kappa
Clonality	Monoclonal
Uniprot ID	Q969J5
Hazard Information	This IL22RA2 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human adrenal gland...



IHC staining of FFPE human adrenal gland...