

---

## Product Datasheet

**IRF4 Antibody / Interferon regulatory factor 4 /  
MUM-1 (orb1823111)**

**Description**

Interferon regulatory factor 4, also called Multiple

ge...

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested Applications**

IHC-P

**Immunogen**

A recombinant partial protein sequence from the C-terminal region of the human protein was used as the immunogen for the MUM1 antibody.

**Storage**

Aliquot the MUM1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Formula**

1 mg/ml in 1X PBS; BSA free, sodium azide free

**Isotype**

Rabbit IgG, kappa

**Clonality**

Recombinant

**Uniprot ID**
[Q15306](#)
**Hazard Information**


This MUM1 antibody is available for research use only.

**Dilution Range**


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

**Expiration Date**


12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...



Negative control: IHC testing of FFPE hu...