
Product Datasheet

**IRF4 Antibody / Interferon regulatory factor 4 /
MUM-1 (orb1823113)**

Description

Interferon regulatory factor 4, also called Multiple

ge...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

Immunogen

A recombinant partial protein sequence from the C-terminal region of the human protein was used as the immunogen for the MUM1 antibody.

Storage

Aliquot the MUM1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Uniprot ID

[Q15306](#)

Hazard Information


This MUM1 antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date


12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



IHC staining of FFPE human tonsil tissue...



Negative control: IHC testing of FFPE hu...