

Product Datasheet

**IRF4 Antibody / Interferon regulatory factor 4 /
MUM-1 (orb1823114)**

Description

Interferon regulatory factor 4 (IRF4), also called



nucl...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P


IHC staining of FFPE human tonsil tissue...

Immunogen

A recombinant partial protein sequence from the C-terminal region of the human protein was used as the immunogen for the IRF4 antibody.

Storage

Aliquot the IRF4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.



Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Rabbit IgG, kappa

Clonality

Recombinant

Uniprot ID

Q15306

Hazard Information

This IRF4 antibody is available for research use only.


Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

IHC staining of FFPE human tonsil tissue...

Expiration Date

12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...