

Product Datasheet

**IRF4 Antibody / Interferon regulatory factor 4 /
MUM-1 (orb1823115)**

Description

Interferon regulatory factor 4 (IRF4), also called



nucl...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P


IHC staining of FFPE human tonsil tissue...

Immunogen

A recombinant partial protein sequence from the C-terminal region of the human protein was used as the immunogen for the IRF4 antibody.

Storage

Aliquot the IRF4 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.



Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Rabbit IgG, kappa

IHC staining of FFPE human tonsil tissue...

Clonality

Recombinant

Uniprot ID


Q15306

Hazard Information

This IRF4 antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT



Expiration Date

12 months from date of receipt.

IHC staining of FFPE human tonsil tissue...