


Product Datasheet

Oct4 Antibody / POU5F1 (orb1823151)

Description	OCT4 (Octamer-binding protein 4), also called POU5F1, is a member of the POU family of transcription...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of the human protein was used as the immunogen for the POU5F1 antibody.
Storage	Aliquot the POU5F1 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	Q01860
Hazard Information	This POU5F1 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human seminoma tiss...