



## Product Datasheet

### IL-3 Antibody (orb1823333)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Interleukin-3, or IL-3, is a pleiotropic cytokine that is primarily secreted by activated T lymphocy... |
| <b>Species/Host</b>        | Mouse                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                            |
| <b>Tested Applications</b> | IHC-P                                                                                                   |
| <b>Immunogen</b>           | A recombinant fragment of human Interleukin 3 protein was used as the immunogen for the IL-3 antibody.  |
| <b>Storage</b>             | Aliquot the IL-3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.       |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Formula</b>             | 1 mg/ml in 1X PBS; BSA free, sodium azide free                                                          |
| <b>Isotype</b>             | Mouse IgG1, kappa                                                                                       |
| <b>Clonality</b>           | Monoclonal                                                                                              |
| <b>Uniprot ID</b>          | <b>P08700</b>                                                                                           |
| <b>Hazard Information</b>  | This IL-3 antibody is available for research use only.                                                  |
| <b>Dilution Range</b>      | Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                         |



Analysis of a  
HuProt (TM)  
microarray  
con...



IHC staining  
of FFPE  
human  
adrenal  
gland...