



Product Datasheet

IL-3 Antibody (orb1823335)

Description	Interleukin-3, or IL-3, is a pleiotropic cytokine that is primarily secreted by activated T lymphocy...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant fragment of human Interleukin 3 protein was used as the immunogen for the IL-3 antibody.
Storage	Aliquot the IL-3 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P08700
Hazard Information	This IL-3 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



Analysis of a HuProt (TM) microarray con...



IHC staining of FFPE human adrenal gland...