



Product Datasheet

IL15 Antibody (orb1823404)

Description	Interleukin-15 (IL-15), also designated IL-T, is a cloned cytokine which shares several biological a...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 1-152) from the human protein was used as the immunogen for the IL15 antibody.
Storage	Aliquot the IL15 antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Recombinant
Uniprot ID	P40933
Hazard Information	This IL15 antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human adrenal gland...



SDS-PAGE analysis of purified, BSA-free ...