
Product Datasheet

**LEF1 Antibody / Lymphoid enhancer-binding factor
1 / TCF1 alpha (orb1823438)**

Description

Lymphoid Enhancing Factor 1 (LEF1), also called



t...

Species/Host

Mouse

Reactivity


Human

Conjugation

Unconjugated

Tested Applications

IHC-P



IHC staining of FFPE human tonsil tissue...

Immunogen

A recombinant fragment of human LEF1 protein was used as the immunogen for the TCF1 alpha antibody.

Storage

Aliquot the TCF1 alpha antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

1 mg/ml in 1X PBS; BSA free, sodium azide free

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Uniprot ID

[Q9UJU2](#)

Hazard Information


This TCF1 alpha antibody is available for research use only.

Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



Analysis of a HuProt (TM) microarray con...