
Product Datasheet

**LEF1 Antibody / Lymphoid enhancer-binding factor
1 / TCF1 alpha (orb1823439)**

Description

Lymphoid Enhancing Factor 1 (LEF1), also called



t...

Species/Host

Mouse

Reactivity

Human

Conjugation


Unconjugated

Tested Applications

IHC-P

Immunogen

A recombinant fragment of human LEF1 protein was used as the immunogen for the TCF1 alpha antibody.



IHC staining of FFPE human tonsil tissue...

Storage

Aliquot the TCF1 alpha antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.

Note

For research use only

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG1, kappa

Clonality

Monoclonal

Uniprot ID

[Q9UJU2](#)

Hazard Information


This TCF1 alpha antibody is available for research use only.

Dilution Range


Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT

Expiration Date

12 months from date of receipt.



IHC staining of FFPE human tonsil tissue...



Analysis of a HuProt (TM) microarray con...