



Product Datasheet

c-Kit Antibody / CD117 / SCFR (orb1823452)

Description	This mAb recognizes a protein of 145kDa, identified as CD117/p145kit. It is found on a wide variety ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A recombinant partial protein sequence (within amino acids 50-250) from the human protein was used as the immunogen for the c-Kit antibody.
Storage	Aliquot the c-Kit antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Rabbit IgG, kappa
Clonality	Recombinant
Uniprot ID	P10721
Hazard Information	This c-Kit antibody is available for research use only.
Dilution Range	ELISA: 2-4ug/ml (Order BSA-free format for coating),Immunohistochemistry (FFPE): 1-2ug/ml for 30 min at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human GIST with c-K...



SDS-PAGE analysis of purified, BSA-free ...